

## *Anti-burn Hydrogel Blanket Kit*



### **Content:**

- N.01 Anti-burn hydrogel blanket of cm 120 x 160
- N.01 Isothermal silver/gold blanket of cm.210 x 160

### **Medical Devices Characteristics:**

- Class CE: II B “ Medical device would be used on use of I, II and III also in presence of injured dermis.
- Features of Burnshield’s medical devices

### **Properties and Advantages:**

- It lowers the temperature of the burned area and rehydrates optimally.
- Helps prevent infections.
- It facilitates the regrowth of the skin.
- The use of Burnshield until full re-epithelialization reduces the possibility of hypertrophic scars

### **Covered Support:**

- Cotton.

### **Technical Specifications:**

The blanket is packaged in an extremely durable and washable synthetic material, with the shape of a double zip closure cylinder. The dimensions of the cover are 210x160, it's in silver and gold colors.

Inside it, as a second protection against any type of impact and puncture, there is a plastic blade of thickness of 1 mm.

The case, equipped with a strap for the shoulder strap, is easily transportable as to leave your hands and all movements free including that of climbing portable stairs and is of red color (V.F. type).

The instructions printed on it are clear, simple and indelible, they are also easy to use by non-practitioners, thanks to the use of figurative instructions that provide any eventuality.

The internal packaging system guarantees rapid use and absolute sterility of the contents. The sterility of the product is guaranteed by an externally made p.v.c. casing and internally in aluminum material.



## Hydrogel Characteristics:

- Idrogel: 96% Water, Melaleuca Alternifolia, Emulsifiers.
- Sterility: Gamma ray irradiation
- PH Gel: 5,5 - 7,0.
- Idrogel Aspect: liquid, white - milky.
- Perfume Idrogel: slightly pharmaceutical.
- Idrogel Biodegradability: 100%.
- Expansion and Storage: 5 years from date of production engraved on each package.  
Store from 15° to 24°.
- Idrogel Incompatibility: none.

## Warnings:

Do not use if the package is damaged or after the expiration date.  
Consult a doctor in case of severe burns.  
In case of intolerance wash thoroughly and consult a doctor.

**Code 50.0701.002**